

- Maximum Brightness 110,000 lx. (LLMX566D model, Battery powered @ 40 cm)
- Single LED X566 have 61,000 lx. (Distance @ 40 cm away)
- Magnetic extension cable: Easy to attach & detach, comfortable to store
- Comes with light-weight handy carrying case
- Battery running time maximum 24 hours (at the lowest level)

Carrying Case



Additional battery
(INCLUDED)



Remote
Controller



FEATURES

- Super bright Dual LED headlight even it's battery powered with ONE control unit
- Compact & light weight control unit, LED module
- Improved & compact size battery pack & Secure assembling system
- LED heat conduction effect: Multi link & attachment is made with aluminum
- LED module weight 12.3 g (LLMX582)
- More soft & durable LED wires