

Control Unit



Carrying Case



66 mm Dual



82 mm Dual



Additional battery
(INCLUDED)



SPECIFICATION

Lens Module Option	Brightness (lx) Min-Max	Operating Time (hr)	Light Field (at 400mm)
LLMX682	3,600 - 41,000	27 - 6 hr	Ø82 mm
LLMX674	4,300 - 47,000	27 - 6 hr	Ø74 mm
LLMX666	5,400 - 61,000	27 - 6 hr	Ø66 mm
LLMX682D	7,100 - 75,000	14 - 2.5 hr	Ø82 mm
LLMX666D	9,500 - 110,000	14 - 2.5 hr	Ø66 mm

FEATURES

- Super bright & lightweight wireless LED Headlight System – Easy to wear
- Perfect compatibility with Headband type loupes
- Eye protection function – Automatic light-off to protect nearby staff's eyes
- Compact battery pack with easy attachment & detachment
- Up to 27 hours of operating time (With single LED lens module)