


"A REAL ADDED VALUE TO
YOUR INTERVENTIONAL
ARMENTARIUM."



The PTCA Balloon solution for challenging lesions

- Low profile, hydrophilic coated balloon catheter to access even more challenging tight lesions
- 21 atm rated burst pressure (RBP)
- 1 single middle marker

a) Lesion entry profile: 0.016"

b) Crossing profile: 0.0205"

c) Hydrophilic coating on balloon and distal shaft

SIS  MEDICAL

Swiss | Interventional | Systems

NIC 1.1[®] hydro

Ordering Information
NIC 1.1[®] hydro CTO-Balloons (RX)

Article Number	Balloon Diameter	Length	RBP	Marker	Min. Guiding
110-006-134	11 mm	6 mm	21 atm	1 single middle	5F
110-010-134	11 mm	10 mm	21 atm	1 single middle	5F
110-015-134	11 mm	15 mm	21 atm	1 single middle	5F
110-020-134	11 mm	20 mm	21 atm	1 single middle	5F





Swiss | Interventional | Systems

Hungerbühlstrasse 12a | CH-8500 Frauenfeld
Tel. +41 (0)52 245 09 90 | Fax +41 (0)52 245 09 91
info@sis-medical.com | www.sis-medical.com

Copyright © by SIS Medical AG
Currently not available for sale in USA.
Subject to alteration due to technical advances.
All rights reserved.
All products are developed and produced by
SIS Medical AG.