ILLUCO

Beyond the vision

LucoCAM LC

ILC-1 Surgical Camera

SEE MORE POSSIBILITIES


4K CLARITY
FOR A CLEARER TOMORROW




- 4K UHD Recording
- Interchangable Lens Options
  2.0x / 3.0x / 5.0x
- Lightweight & Compact
- Easy Connection
  USB Type-C
- Bright & Clear Imaging

ILLUCO
Beyond the vision


Smart Recording Made Simple




Easy Recording & Capture




Smart Bookmarking



Freeze & Digital Zoom



Quick Video Trimming



Intuitive Library & Playback




Flexible Control & Versatile Accessories


Flexible

Mounting Options

Compatible with both headband and loupe setups.
Also compatible with gooseneck
and third-party mounts.



Loupe Mount
(Bolt-on / Magnetic)



Headband Mount
(Soft / Hard)

Technical Specifications

Resolution	UHD(4K,2160p), FHD (1080p), HD (720p)
Lens Magnification	2.0X / 3.0X , 5.0X (Option)
Focus Mode	Manual Focus
Working distance	40-100 cm
Audio Input	Digital Microphone
Output & Power	High Speed USB 2.0 (DC 5V, 270mA)
Supported OS	Windows, Mac, Android, iPadOS
Dimensions / Weight	23 × 23× 38.7 mm / 35g

Software and Manual Download

www.illuco.co.kr
Support > Download

ILLUCO
Beyond the vision
www.illuco.co.kr T: +82-31-429-8825 E: info@illuco.co.kr